

Elos Accurate[®] IO Scan Body

For an Open Digital Workflow



**One screwdriver
regardless of implant
system**

**Color coded according to
implant manufacturer and
platform**

**PEEK scanning part for best
scanning result**



Our products are purchased through our web:
<https://emeralddental.com/>


1-888-251-7171

Elos Accurate® Analog for Printed Models


To ensure a reliable model

The new Elos Accurate Analog for Printed Models is specially made for desktop 3D printers. All it's specific features are specified here:


Available for the major implant platforms and color-coded according to implant manufacturer



3 strong locking areas ensure the correct position




Inside the model:
Undercut in print socket



Elos Accurate Analog for Printed Models (PMA)

- Precise height installation
- Secure press-fit in the model
- One insertion position — always correct!
- Automatically centered
- Complete with installation tools — one set for all platforms



Our popular analog Elos Accurate Model Analog (MA) is still available

- Designed to fit milled or plaster models
- Audible click sound confirms correct installation in the model
- Single-use for predictable outcome
- Can be used for plaster, milled and printed models.

Our products are purchased through our web:

<https://emeralddental.com/>

1-888-251-7171